


DIGNISHIELD[®]

STOOL MANAGEMENT SYSTEM

Stool Management System Maintenance Algorithm

This document provides a quick reference of maintenance activities related to the DIGNISHIELD[®] Stool Management System (SMS). For further information, please consult the product instructions for use.




Bard and DigniShield are trademarks and/or registered trademark of C. R. Bard, Inc. Please consult product labels and inserts for any hazards, warnings, cautions and directions for use. ©2013 C. R. Bard, Inc. All Rights Reserved. 1307-37 R10/13 THP P10/13 XM

Directions for use

ATTENTION: PROPER INSERTION IS IMPORTANT

This card provides a quick reference for the preparation and insertion of the DIGNISHIELD® SMS Cuff. Consult product IFU for further information on the use of this device.

- 1 DEFLATE**
Remove all air from cuff using syringe.




- 2 FLATTEN/FOLD**
Flatten cuff, fold the upper right corner down to the bottom of the cuff in a 45° angle. This creates a leading edge for easy insertion.




- 3 LUBRICATE**
Apply lubrication.




- 4 INSERT**
Insert cuff into rectal vault.



- 5 INFLATE (45ml water)**
Inflate cuff with syringe and seat.



- 6 TEST**
Irrigate to verify seating.



The Bard® DigniShield® Stool Management System (SMS) is a single-use device intended for fecal management by diverting and collecting liquid or semi-liquid stool to minimize skin contact in adult, bedridden patients.

Contraindications: The device should not be used for more than 29 consecutive days, on patients with certain medical conditions including rectal or anal injury or abnormalities, on patients with rectal mucosa impairment or on patients who have had lower large bowel or rectal surgery within the last year. Do not use on patients with indwelling rectal or anal devices, delivery mechanisms, or enemas in place.

Warnings and Adverse Events: As with the use of any rectal device, adverse events can occur including: leakage of stool, loss of anal sphincter muscle tone, pressure necrosis, infection, bowel obstruction, and perforation of the bowel. Changes to the patient including: rectal bleeding indicating possible pressure necrosis, abdominal distention, cuff migration, and rectal pain should be investigated.

Please consult package insert for more detailed safety information and instructions for use.

BARD

8195 Industrial Blvd.
Covington, GA 30014
1.800.526.4455
www.bardmedical.com

If you have questions or require additional information, please contact your local Bard Sales Representative or call 1-800-526-4455.

Bard and DigniShield are trademarks and/or registered trademarks of C. R. Bard, Inc.

©2013 C. R. Bard, Inc. All Rights Reserved.