

EXTENDED CARE UROLOGY PRODUCTS

EXTENDED CARE

UROLOGY PRODUCTS



UROLOGY

BAIRD | MEDICAL

BARDIA® URETHRAL CATHETERS

Red rubber all-purpose urethral catheters:

- Radiopaque red rubber catheter with a round, hollow tip
- Can be used as a Robinson or a Nelaton catheter
- Two opposing drainage eyes
- 16" long
- Single-use, sterile

802408-22	8-22 Fr.	100/case
-----------	----------	----------



VINYL CATHETERS

BARD® CLEAN-CATH® Intermittent catheter:

- Finger grips for ease of use
- Flexible funnel
- Latex-Free
- Individually packaged in easy to open peel-strip package
- Smooth round tip for easy insertion
- Olive coude tip helpful for insertion difficulties
- Single-use, sterile

6 in. long	420706-18	6-18 Fr.	50/case
10 in. long	421108-10	8-10 Fr.	50/case
16 in. long	421708-18	8-18 Fr.	50/case
Coude Tip	423710-18	10-18 Fr.	30/case

without funnel

6 in. long	420614	14 Fr.	50/case
10 in. long	421008-10	8-10 Fr.	50/case
16 in. long	421610-14	10-14 Fr.	50/case



BARDIA® URETHRAL TRAYS

Complete urethral catheterization tray:

- Waterproof underpad, fenestrated drape, 3 BZK swabsticks, 5g lubricating jelly, 2 latex-free exam gloves, specimen container with label
- Peel-top tray
- Single-use, sterile

802115	With 15 Fr. red rubber urethral catheter	20/case
802214	With 14 Fr. plastic urethral catheter	20/case
802100	Without catheter	20/case



TOUCHLESS® PLUS UNISEX INTERMITTENT CATHETER KIT

- Closed system kit allows for ease of drainage, simplifies disposal and reduces the risk of cross-contamination
- Pre-lubricated catheter with finger guide
- 1100cc collection chamber with built-in sample/drainage port
- 2 latex-free exam gloves, povidone iodine swabs and underpad
- Single use, sterile

Vinyl- Latex-Free

4A5108, 10	8, 10 Fr	50/case
4A5142, 44	12, 14 Fr.	50/case
No accessories included		
4A6110, 42, 44	10, 12, 14 Fr	50/case

Red rubber also available

4A5042, 44	12, 14 Fr.	50/case
4A7044 (coude)	14 Fr.	50/case
No accessories included		
4A6044	14 Fr.	50/case

BARDIA® FOLEY INSERTION TRAYS

Foley catheter prepping components only:

- Waterproof underpad, fenestrated drape, 5g lubricating jelly, 2 latex-free exam gloves, specimen container with label and/or graduated collection container
- For 5cc or 30cc balloon Foley catheters
- 3 povidone iodine or 3 BZK swabsticks
- CSR-wrap package or Peel-top package
- Single-use, sterile

Six Versions Available - See Back Page




DRAINAGE BAGS

Drainage Bags with swivel hanger/ flexible hook and loop cloth strap which facilitate attachment to bed frame or wheelchair:

- 2000ml capacity with anti-reflux chamber
- Urine sampling area
- Easy-to-drain outlet tube with clamp and tube holder
- Plastic bag, single-use, sterile fluid path
- Latex-Free

802001	Urinary Drainage Bag swivel hanger with flexible hook	20/case
802002	Urinary Drainage Bag with flexible hook and loop cloth strap	20/case



COATED LATEX FOLEY CATHETERS

Silicone-elastomer coated, 2-way, latex Foley catheters:

- Silicone-elastomer coating provides protection against urethral irritation
- Available with either a 5cc or a 30cc balloon
- Peel-pouch
- Single-use, sterile

123512A-30A	12-30 Fr., 5cc	12/case
123612A-30A	12-30 Fr., 30cc	12/case

ALL-SILICONE FOLEY CATHETERS

100% silicone, 2-way, Foley catheters:

- Available with either a 5cc or a 30cc balloon
- Peel-pouch
- Single-use, sterile
- Latex-Free

806514-24	14-24 Fr., 5cc	12/case
806316-24	16-24 Fr., 30cc	12/case



PISTON SYRINGE

750379	BARD® Enteral 60cc Piston syringe in polyethylene bag, non-sterile	30/case
--------	---	---------



BARD® DISPOZ-A-BAG® LEG BAG AND ACCESSORIES

- The industry standard for quality, reliability and dependability
- Made of heavy duty vinyl for long wear time, odor containment, and leak resistance
- FLIP-FLO™ drainage valve opens and closes with a simple flip
- Reusable, sterile



Leg Bag 4-Pack

150419, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (Bonus pack 4 leg bags, 1 pair fabric leg straps)	12/case
------------	--	---------

Leg Bag Systems

150819, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair fabric leg straps, 18" PVC extension tubing)	50/case
150919, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (18" PVC extension tubing, 4 1/2" drainage tube, no straps)	50/case
151919, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (18" PVC extension tubing, 4 1/2" drainage tube, 1 pair latex straps)	50/case
650819, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair fabric leg straps, 18" latex extension tubing)	50/case
651919, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (18" latex extension tubing, 4 1/2" drainage tube, no straps)	50/case

Leg Bag with Fabric Straps

150719, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair fabric straps)	50/case
------------	---	---------

Leg Bag with Latex Straps

150106, 107, 108	9, 19, 32 oz. Leg Bag with rubber cap valve (1 pair latex straps)	48/case
150319, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair latex straps)	50/case
4A4173, 74	19, 32 oz. Leg Bag with drain tube and clamp (1 pair latex straps)	24/case

Leg Bag with No Straps

150619, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (no straps)	50/case
------------	--	---------

Straps

150707	Latex leg straps - 22" long	50 pair/ case
150507	Deluxe fabric leg straps - 24" long, 3/4" wide	24 pair/ case
162110 (Small) 162210 (Medium) 162310 (Large)	2" wide leg straps with fabric backing and Velcro™ closure (9" - 13" long, Small) (13" - 20" long, Medium) (20" - 27" long, Large)	6 pair/ case

Tubing

00538	8 1/2" tubing and adapter to connect leg bag to male external or Foley catheters	50/case
4A4194	18" PVC Extension tubing with connector, Sterile	50/case
6A4194	18" Latex Extension tubing with connector, Sterile	50/case

Leg Bag Holders

150111 (Small) 150121 (Medium) 150131 (Large)	Leg Bag Holders (8" - 14", Small) (14" - 20", Medium) (16" - 26", Large)	12/case
---	---	---------

Bard EXTENDED CARE Urology Products

Reorder Code	Foley Catheter	Drain Bag (Cat. No. 802001)	Prepping Components				Package		Quantity per case	HCPCS* Code
			PVI Swabs	BZK Swabs	10cc Syringe	30cc Syringe	CSR Wrap	Peel-Top		
CLOSED SYSTEM FOLEY TRAYS										
802016**	16 Fr. 5cc	✓	✓		✓		✓	10	A4314	
802018**	18 Fr. 5cc	✓	✓		✓		✓	10	A4314	
ADD-A-FOLEY TRAYS										
802015		✓	✓		✓		✓	10	A4354	
802035		✓	✓		✓		✓	10	A4354	
FOLEY INSERTION TRAYS										
802010			✓		✓		✓	20	A4310	
802030			✓				✓	20	A4310	
802110				✓	✓		✓	20	A4310	
802130				✓		✓	✓	20	A4310	
802011			✓		✓		✓	20	A4310	
802031			✓		✓		✓	20	A4310	
FOLEY CATHETERS										
1235XXA	5cc balloon, silicone elastomer coated latex Foley catheter (12-30 Fr.)***							12	A4338	
1236XXA	30cc balloon, silicone elastomer coated latex Foley catheter (12-30 Fr.)***							12	A4338	
8065XX	5cc balloon, all-silicone Foley catheter (14-24 Fr.)***							12	A4344	
8063XX	30cc balloon, all-silicone Foley catheter (16-24 Fr.)***							12	A4344	
DRAINAGE BAGS										
802001	2 liter urine drainage bag w/ anti-reflux chamber and swivel hanger w/ flexible hook							20	A4357	
802002	2 liter urine drainage bag w/ anti-reflux chamber and flexible hook and loop cloth strap							20	A4357	
URETHRAL CATHETERS & TRAYS										
802115	15 Fr. red rubber urethral catheter, lube, BZK swabs, gloves, drapes, specimen jar							20	A4353	
802100	Same as above without catheter							20		
802214	14 Fr. vinyl urethral catheter, lube, BZK							50	A4353	
8024XX	Red Rubber urethral catheter (8-22 Fr.)***							50	A4351	
4207XX	Vinyl catheter with funnel - 6 inches (6-18 Fr.)***							50	A4351	
4211XX	Vinyl catheter with funnel - 10 inches (8-10 Fr.)***							50	A4351	
4217XX	Vinyl catheter with funnel - 16 inches (8-18 Fr.)***							50	A4351	
4237XX	Olive coude tip vinyl catheter w/ funnel - 16 inches (10-18 Fr.)***							30	A4352	
4A5142	Closed system kit with vinyl catheter (12 Fr.)							50	A4353	
4A5144	Closed system kit with vinyl catheter (14 Fr.)							50	A4353	
4A5042	Closed system kit w/ Red Rubber catheter (12 Fr.)							50	A4353	
4A5044	Closed system kit with Red Rubber catheter (14 Fr.)							50	A4353	
4A7044	Closed system kit with Red Rubber coude catheter (14 Fr.)							50	A4353	
EXTENSION TUBING										
4A4194	18" Latex-Free Extension Tubing, sterile							50		
6A4194	18" Latex Extension Tubing, sterile							50		

✓ - Indicates that the item is included in tray. Prepping components include: Swabsticks as indicated in the table, pre-filled inflation syringe as indicated in the table, lube jelly, exam gloves, fenestrated drape, waterproof underpad, and a graduated collection container (except where noted).

**In addition to prepping components, a specimen container, a cap and label are also included.

***Where XX is equal to the French size (even Fr. sizes only)

ATTENTION: All products featured in this brochure are available for ordering, however not all products are listed in the tables.

For technical or clinical information on any product, call Bard's Medical Services and Support (MSS) staff, toll-free at 1-800-227-3357. Physician-to-Physician consultation is also available for patient-specific inquiries. For a complete listing of all Extended Care products, contact Bard Medical Division or your local sales representative.

*NOTE: Reimbursement policies vary widely from insurer and reimbursement policies of the same insurer may vary in different sections of the U.S. As reimbursement policies are subject to change, the information contained herein may not be accurate at the time of use. Prior to submission of a claim for reimbursement, the user should contact the insurer (i.e., Medicare, Medicaid, or private payer) to verify applicable codes and reimbursement levels. C. R. Bard, Inc. does not guarantee that the purchases of any similar item is reimbursable in whole or in part, by any public or private health care insurer. Source for HCPCS codes: <http://www.hcfa.gov/stats/pfiles.htm>

BARD | MEDICAL

C. R. Bard, Inc.
Covington, GA 30014
1-800-526-4455
www.bardmedical.com

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use. Bard, Bardia, Flip-Fllo, Touchless, Clean-Cath and Dispoz-a-Bag are trademarks and/or registered trademarks of C. R. Bard, Inc.

©2011 C. R. Bard, Inc. All Rights Reserved. 1202-17 R03/12 THP P03/12 5M