


INTERMITTENT CATHETER SYSTEM

TOUCHLESS[®] PLUS

INTERMITTENT CATHETER SYSTEM



Designed to help reduce the risk
of urinary tract infections.

UROLOGY

BAIRD | MEDICAL


Safe, Sterile and Simple

The BARD® TOUCHLESS® Plus Intermittent Catheter System

The BARD® TOUCHLESS® Plus Intermittent Catheter system has what is known as a “closed system” designed to help reduce the risk of urinary tract infections. Among the features that you’ll find helpful are:

- 1 Insertion tip** correctly positions the catheter in the urethra to help reduce UTIs.
- 2 Catheter control guide** allows control of catheter insertion with the touch of a finger.
- 3 Tapered neck design** decreases cuts and scrapes associated with sharp corners.
- 4 Catheter cap** has pull ring for easy removal.
- 5 Pre-lubricated catheter** reduces discomfort during insertion.
- 6 Large 1100cc** collection bag is durable, virtually eliminating the risk of leaks.

Complete instructions for use included in product packaging.



Easy to Use

Easy to insert, easy to use, easy to remove - the features needed to enhance convenience and usage.

- Catheter cap has pull ring for easy removal
- Pre-lubricated catheter reduces discomfort during insertion.
- Large 1100cc collection bag is durable, virtually eliminating the risk of leaks.



TOUCHLESS® PLUS

INTERMITTENT CATHETER SYSTEM

Your System, Your Choice

- Choice of red rubber or vinyl catheters based upon your preference and needs.
- Variety of catheter sizes, and available with or without accessories.



Product Number	Description	French Size	Latex-Free
Accessories included - 3 PVI swabs, 2 latex-free gloves, underpad			
4A5108	Vinyl	8 Fr.	Yes
4A5110	Vinyl	10 Fr.	Yes
4A5142	Vinyl	12 Fr.	Yes
4A5144	Vinyl	14 Fr.	Yes
NEW 4A5146	Vinyl	16 Fr.	Yes
4A7114	Vinyl Coude	14 Fr.	Yes
4A7116	Vinyl Coude	16 Fr.	Yes
4A5042	Red Rubber	12 Fr.	No
4A5044	Red Rubber	14 Fr.	No
4A7044	Red Rubber Coude	14 Fr.	No

For more information, please contact us at 1-800-526-4455 or visit www.bardmedical.com.

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use.

BARD | MEDICAL

ESSENTIAL. EFFECTIVE. EFFICIENT.

C. R. Bard, Inc.
Covington, GA 30014
www.bardmedical.com

Bard and Touchless are registered trademarks of C. R. Bard, Inc.
©2012 C. R. Bard, Inc. All Rights Reserved. Printed in U.S.A.
1207-31 R07/12 THP P07/12 5M