

200SMS REF SMS002
BARD | MEDICAL

DIGNISHIELD
STOOL MANAGEMENT SYSTEM



BARD | MEDICAL
REF SMS002

STOOL MANAGEMENT SYSTEM
DIGNISHIELD

DIGNISHIELD™

STOOL MANAGEMENT SYSTEM

Additional DIGNISHIELD™ Stool Management System Collection Bags available in 10 count cartons
Order number SMS2B1L

REF SMS002

The device has no components made of natural rubber latex.



Bard and Dignishield are trademarks and/or registered trademarks of C. R. Bard, Inc.

Manufactured by:
BARD
C. R. Bard, Inc.
Covington, GA 30014, U.S.A.
1.800.526.4455
www.bardmedical.com
Bard Medical Division
Assembled in Mexico
PK7634884 12/2011

Contents:



ASSESS PATIENT SELECTION AND SMS CONTRAINDICATIONS BEFORE PROCEEDING

Patient Selection Parameters for SMS Use: Check All That Apply

- Rectal Tone: Confirmed by Digital Rectal Exam
- Fecal Impaction: Absent
- Stool Frequency: 3-6 Stools in 24 hrs
- Stool Consistency: Semi-Liquid to Liquid
- Bristol Stool Form Scale: Type 6-7 (diarrhea)
- Stool Culture: Positive for *C. difficile*
- Braden Scale Score: ≤ 14
- IAD* Risk Level: High Risk to Severe IAD

*Incontinence Associated Dermatitis

CONTRAINDICATIONS for SMS Use: If Any Box Checked DO NOT PROCEED

- Pediatric patient
- Previous uninterrupted use of SMS for the past 29 days
- Known sensitivity or allergy to any system components
- Lower large bowel or rectal surgery within last year
- Patients with any rectal or anal injury, stricture, stenosis, tumor, severe hemorrhoids or fecal impactions
- Suspected or confirmed rectal mucosal impairment
- Indwelling rectal or anal device, delivery mechanism or enemas in place
- Rectal bleeding precautions

BARD | MEDICAL
REF SMS002

DIGNISHIELD
STOOL MANAGEMENT SYSTEM

BARD | MEDICAL

DIGNISHIELD™
STOOL MANAGEMENT SYSTEM

BARD | MEDICAL
REF SMS002